



Product Datasheet

IKBKB antibody (orb400952)


Description	IKBKB antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Inhibitor of nuclear factor kappa-B kinase subunit beta protein (395-588AA)
Target	IKBKB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O14920
Dilution Range	WB:1:500-5000, IHC-P:1:500-1000, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 whole cell...



Immunohistochemical staining of human br...



Immunofluorescence analysis of A549 cell...