



Product Datasheet

EIF3F antibody (orb40104)

Description	EIF3F antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP
Immunogen	Recombinant Human Eukaryotic translation initiation factor 3 subunit F protein (88-357AA)
Target	EIF3F
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	000303
Dilution Range	IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human brain tissue.



Immunohistochemical staining of human adrenal gland tissue.