


Product Datasheet

Pectate lyase 1 antibody (orb401100)

Description	Pectate lyase 1 antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Hesperocyparis arizonica Pectate lyase 1 protein (22-367AA)
Target	Pectate lyase 1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9SCG9
Dilution Range	WB:1:500-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of Recombinant pro...