

Product Datasheet

Pectate lyase 1 antibody (FITC) (orb401102)

Description	Pectate lyase 1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	FITC
Immunogen	Recombinant Hesperocyparis arizonica Pectate lyase 1 protein (22-367AA)
Target	Pectate lyase 1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9SCG9
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.