



Product Datasheet

TGFB111 antibody (orb401268)

Description	TGFB111 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Transforming growth factor beta-1-induced transcript 1 protein (138-200AA)
Target	TGFB111
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O43294
Dilution Range	IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293 whole ce...



Immunohistochemical staining of human ce...