

Product Datasheet

CSF2RB antibody (Biotin) (orb401455)

Description	CSF2RB antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Cytokine receptor common subunit beta protein (492-710AA)
Target	CSF2RB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P32927
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.