

Product Datasheet

Pollen allergen Phl p 5a antibody (orb401688)

Description	Pollen allergen Phl p 5a antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Phleum pratense Pollen allergen Phl p 5a protein (1-286AA)
Target	Pollen allergen Phl p 5a
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q40962
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.