

Product Datasheet

Pollen allergen Phl p 5a antibody (HRP) (orb401689)

Description	Pollen allergen Phl p 5a antibody (HRP)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Phleum pratense Pollen allergen Phl p 5a protein (1-286AA)
Target	Pollen allergen Phl p 5a
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q40962
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.