



Product Datasheet

WWOX antibody (orb401748)


Description	WWOX antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	Recombinant Human WW domain-containing oxidoreductase protein (1-180AA)
Target	WWOX
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NZC7
Dilution Range	WB:1:500-5000, IHC-P:1:200-500, IF:1:50-200, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of MCF-7 whole cel...



Immunohistochemical staining of human te...



Immunohistochemical staining of human pl...