

Product Datasheet

JD5037 (orb401799)

Description	JD5037 is an inverse agonist of peripheral cannabinoid 1 (CB1) receptors (K _i = 0.35 nM).
Note	For research use only
Formula	C ₂₇ H ₂₇ Cl ₂ N ₅ O ₃ S
MW	572.506
CAS Number	[1392116-14-1]
SMILES	<chem>C1C([H])=C([H])C(=C([H])C=1[H])C1[C@@]([H])(C2C([H])=C([H])C([H])=C([H])C=2[H])C([H])([H])N/C/N([H])S(C2C([H])=C([H])C(=C([H])C=2[H])Cl)(=O)=O)=N[C@]([H])(C(N([H])[H])=O)C([H])(C([H])([H])[H])C([H])([H])[H])N=1</chem>
Expiration Date	12 months from date of receipt.