

Product Datasheet

GNE-3511 (orb401810)

Description	GNE-3511 is a potent and selective zipper kinase (DLK, MAP3K12) inhibitors.
Note	For research use only
Formula	$C_{23}H_{26}F_2N_6O$
MW	440.489
CAS Number	[1496581-76-0]
SMILES	<chem>FC1(C([H])([H])C([H])([H])N(C1([H])[H])C1=C([H])C(=C([H])C(N([H])C2C([H])=C(C#N)C([H])=C([H])N=2)=N1)C1([H])C([H])([H])C([H])([H])N(C([H])([H])C1([H])[H])C1([H])C([H])([H])OC1([H])[H])F</chem>
Expiration Date	12 months from date of receipt.