

Product Datasheet

Olodaterol (orb402067)

Description	Olodaterol (BI1744) is a long acting β 2-adrenoceptor agonist with an EC50 of 1.4 ± 0.08 nM.
Note	For research use only
Formula	$C_{21}H_{26}N_2O_5$
MW	386.44
CAS Number	[868049-49-4]
SMILES	<chem>O=C1NC2=CC(O)=CC([C@@H](O)CNC(C)(C)CC3=CC=C(OC)C=C3)=C2OC1</chem>
Expiration Date	12 months from date of receipt.