

## Product Datasheet

### GSK682753A (orb402074)

<b>Description</b>	GSK682753A is a selective and highly potent inverse agonist of the epstein-barr virus-induced recept...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{23}H_{21}Cl_3N_2O_3$
<b>MW</b>	479.78
<b>CAS Number</b>	[1334294-76-6]
<b>SMILES</b>	<chem>O=C1OC2(CCN(C/C=C/C3=CC=C(CI)C=C3)=O)CC2)CN1CC4=CC=C(CI)C(CI)=C4</chem>
<b>Expiration Date</b>	12 months from date of receipt.