

Product Datasheet

CD437 (orb402171)

Description	CD437 is a synthetic retinoid that is an RAR γ -selective agonist. It displays RAR γ -dependent and -i...
Note	For research use only
Formula	C ₂₇ H ₂₆ O ₃
MW	398.49354
CAS Number	[125316-60-1]
SMILES	<chem>O=C(C1=CC=C2C=C(C3=CC=C(O)C(C4(C5)CC6CC5CC(C6)C4)=C3)C=CC2=C1)O</chem>
Expiration Date	12 months from date of receipt.