


## Product Datasheet

### TRIM69 antibody (orb40239)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | TRIM69 antibody  |
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Human, Mouse   |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Tested Applications</b> | ELISA, IHC, WB   |
| <b>Immunogen</b>           | Recombinant Human E3 ubiquitin-protein ligase TRIM69 protein (1-341AA) |
| <b>Target</b>              | TRIM69   |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH7.3.                      |
| <b>Form/Appearance</b>     | Liquid   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.          |
| <b>Note</b>                | For research use only  |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Uniprot ID</b>          | <a href="#">Q86WT6</a>   |
| <b>Dilution Range</b>      | IHC-P:1:20-200   |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |



Western blot analysis of Mouse kidney(la...