


Product Datasheet

Human BMP-9 ELISA Kit (orb402721)

Description	Human BMP-9 ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: S320-R429
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is around 120 pg/ml
Sample Types	bone tissue, cell culture supernatants, serum plasma (heparin, EDTA) and urine.
Uniprot ID	Q9UK05
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human BMP-9 ELISA Kit standard curve