


Product Datasheet

Rat Clusterin ELISA Kit (orb402729)

Description	Rat Clusterin ELISA Kit
Reactivity	Rat
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D22-R226 (β) & S227-E448 (α)
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:20000, concentration in serum and plasma is around 80 ug/ml
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	P05371
Sensitivity	15 pg/ml
Expiration Date	Please enquire.



Rat Clusterin ELISA Kit standard curve