

Product Datasheet

pBLNK antibody (FITC) (orb4028)

Description	Rabbit polyclonal to pBLNK.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthesised phosphopeptide derived from human BLNK around the phosphorylation site of Tyr189 [EN(p-Y)IH]
Target	BLNK
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	65kDa
Uniprot ID	Q8WV28
Dilution Range	IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.