
Product Datasheet

Human Galectin-4 / LGALS4 / L36LBP Protein (orb402834)

Description

Human Galectin-4, His Tag is

(HEK293). It contains AA Met 1 - Ile 323...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Galectin-4

Preservatives

50 mM HEPES, 150 mM NaCl, 5 mM EDTA, 2 mM TCEP, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

N-9×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 38.1 kDa. The protein migrates as 40-44 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

NP_006140.1

Purity

95%

Source

Human Galectin-4, His Tag (orb402834) is expressed from human 293 cells (HEK293). It contains AA Met 1 - Ile 323 (Accession # P56470-1).

MW

38.1 kDa

Uniprot ID[P56470](#)**NCBI**[NP_006140.1](#)**Expiration Date**

6 months from date of receipt.