

---

## Product Datasheet

### Mouse OX40 Ligand / TNFSF4 Protein (orb402996)

**Description**

Mouse OX40 Ligand, His Tag is

(HEK293). It contains AA Gln 49 - Leu 1...

<b>Reactivity</b>	Mouse
<b>Endotoxins</b>	1.0 EU per µg
<b>Conjugation</b>	Unconjugated
<b>Target</b>	OX40 Ligand
<b>Preservatives</b>	PBS, pH7.4
<b>Form/Appearance</b>	Powder
<b>Storage</b>	-20°C
<b>Tag</b>	N-10×His
<b>Note</b>	For research use only
<b>Application notes</b>	This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 18.3 kDa. The protein migrates as 20 kDa and 22-25 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.
<b>Protein Sequence</b>	NP_033478.1
<b>Purity</b>	90%
<b>Source</b>	Mouse OX40 Ligand, His Tag (orb402996) is expressed from human 293 cells (HEK293). It contains AA Gln 49 - Leu 198 (Accession # P43488-1).
<b>MW</b>	18.3 kDa
<b>Uniprot ID</b>	<a href="#">P43488</a>
<b>NCBI</b>	<a href="#">NP_033478.1</a>
<b>Expiration Date</b>	6 months from date of receipt.