

## Product Datasheet

### Bovine Transthyretin ELISA Kit (orb403314)

<b>Description</b>	Bovine Transthyretin(TTR) ELISA kit
<b>Reactivity</b>	Bovine
<b>Range</b>	312.5 ng/mL-5000 ng/mL
<b>Target</b>	TTR
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>O46375</b>
<b>Sensitivity</b>	312.5 ng/mL
<b>Expiration Date</b>	Please enquire.