

Product Datasheet

Bovine TGF beta 1 ELISA Kit (orb403324)

Description	Bovine transforming growth factor β 1, TGF- β 1 ELISA Kit
Reactivity	Bovine
Range	6.25 ng/mL-400 ng/mL
Target	TGFB1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates
Uniprot ID	P18341
Sensitivity	1.56 ng/mL
Expiration Date	Please enquire.