

Product Datasheet

Dog FSH ELISA Kit (orb403535)

Catalog Number	orb403535
Category	Assays and Kits
Description	The Dog FSH ELISA Kit detects the target FSH in canine samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 0.35 mIU/mL-140 mIU/mL. The kit offers a sensitivity of 0.35 mIU/mL.
Target	FSH
Reactivity	Canine
UniProt ID	E2RJZ8
Application notes	Detection Wavelength: 450 nm
Assay Type	Competitive
Assay Time	1-5h
Range	0.35 mIU/mL-140 mIU/mL
Sensitivity	0.35 mIU/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)