

## Product Datasheet

### Goat CGA ELISA Kit (orb403695)

<b>Description</b>	Goat Thyroid Stimulating Hormone, TSH ELISA Kit
<b>Reactivity</b>	Goat
<b>Range</b>	0.5 µIU/mL-20 µIU/mL
<b>Target</b>	CGA
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">Q8WMW8</a>
<b>Sensitivity</b>	0.25 µIU/mL
<b>Expiration Date</b>	Please enquire.