

## Product Datasheet

### Goat Luteinizing Hormone ELISA Kit (orb403709)

<b>Description</b>	Goat Luteinizing Hormone(LH) ELISA KIT
<b>Reactivity</b>	Goat
<b>Range</b>	0.3 mLU/mL-60 mLU/mL
<b>Target</b>	LH
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Sensitivity</b>	0.15 mIU/mL
<b>Expiration Date</b>	Please enquire.