

## Product Datasheet

### Guinea Pig IL1B ELISA Kit (orb403755)

<b>Description</b>	Guinea pig Interleukin-1 beta(IL1B) ELISA kit
<b>Reactivity</b>	Guinea pig
<b>Range</b>	15.6 pg/mL-1000 pg/mL
<b>Target</b>	IL1B
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">Q9WVG1</a>
<b>Sensitivity</b>	3.9 pg/mL
<b>Expiration Date</b>	Please enquire.