

## Product Datasheet

### Guinea Pig MCP1 ELISA Kit (orb403765)

<b>Description</b>	Guinea pig monocyte chemotactic protein 1, MCP1 ELISA Kit
<b>Reactivity</b>	Guinea pig
<b>Range</b>	23.5 pg/mL-1500 pg/mL
<b>Target</b>	CCL2
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">Q08782</a>
<b>Sensitivity</b>	5.8 pg/mL
<b>Expiration Date</b>	Please enquire.