

Product Datasheet

Human YWHAQ ELISA Kit (orb403875)

Description	Human 14-3-3 protein theta(YWHAQ) ELISA kit
Reactivity	Human
Range	0.47 ng/mL-30 ng/mL
Target	YWHAQ
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cerebrospinal fluid (CSF)
Uniprot ID	P27348
Sensitivity	0.11 ng/mL
Expiration Date	Please enquire.