

Product Datasheet

Human BCL2L11 ELISA Kit (orb404018)

Description	Human Bcl-2-like protein 11(BCL2L11) ELISA kit
Reactivity	Human
Range	23.44 pg/mL-1500 pg/mL
Target	BCL2L11
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	O43521
Sensitivity	5.86 pg/mL
Expiration Date	Please enquire.