

## Product Datasheet

### Human CPE ELISA Kit (orb404082)

<b>Description</b>	Human Carboxypeptidase E(CPE) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	39.06 pg/mL-2500 pg/mL
<b>Target</b>	CPE
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>P16870</b>
<b>Sensitivity</b>	9.76 pg/mL
<b>Expiration Date</b>	Please enquire.