

Product Datasheet

Human CEACAM8 ELISA Kit (orb404093)

| | |
|--------------------------|--|
| Description | Human Carcinoembryonic antigen-related cell adhesion molecule 8(CEACAM8) ELISA kit |
| Reactivity | Human |
| Range | 0.625 ng/mL-40 ng/mL |
| Target | CEACAM8 |
| Note | For research use only |
| Application notes | Detection Wavelength: 450 nm |
| Sample Types | serum, plasma, tissue homogenates, cell lysates |
| Uniprot ID | P31997 |
| Sensitivity | 0.156 ng/mL |
| Expiration Date | Please enquire. |