

## Product Datasheet

### Human COL10A1 ELISA Kit (orb404172)

<b>Description</b>	Human Collagen alpha-1(X) chain (COL10A1) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	0.312 ng/mL-20 ng/mL
<b>Target</b>	COL10A1
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates
<b>Uniprot ID</b>	<a href="#">Q03692</a>
<b>Sensitivity</b>	0.078 ng/mL
<b>Expiration Date</b>	Please enquire.