

## Product Datasheet

### Human DEFA6 ELISA Kit (orb404281)

<b>Description</b>	Human Defensin-6(DEFA6) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156 ng/mL-10 ng/mL
<b>Target</b>	DEFA6
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">Q01524</a>
<b>Sensitivity</b>	0.039 ng/mL
<b>Expiration Date</b>	Please enquire.