

## Product Datasheet

### Human DLL4 ELISA Kit (orb404287)

<b>Description</b>	Human Delta-like protein 4(DLL4) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	23.5 pg/mL-1500 pg/mL
<b>Target</b>	DLL4
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>Q9NR61</b>
<b>Sensitivity</b>	5.8 pg/mL
<b>Expiration Date</b>	Please enquire.