

Product Datasheet

Human ADAM2 ELISA Kit (orb404303)

Description	Human Disintegrin and metalloproteinase domain-containing protein 2(ADAM2) ELISA kit
Reactivity	Human
Range	0.312 ng/ml - 20 ng/ml
Target	ADAM2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q99965
Sensitivity	0.078 ng/ml
Expiration Date	Please enquire.