

## Product Datasheet

### Human EDIL3 ELISA Kit (orb404355)

**Description** Human EGF-like repeat and discoidin I-like domain-containing protein 3(EDIL3) ELISA kit

**Reactivity** Human

**Range** 125 pg/mL-8000 pg/mL

**Target** EDIL3

**Note** For research use only

**Application notes** Detection Wavelength: 450 nm

**Sample Types** serum, plasma, tissue homogenates

**Uniprot ID** [O43854](#)

**Sensitivity** 31.25 pg/mL

**Expiration Date** Please enquire.