

## Product Datasheet

### Human EpCAM ELISA Kit (orb404392)

<b>Description</b>	Human Epithelial cell adhesion molecule(EPCAM) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	0.625 ng/ml-40 ng/ml
<b>Target</b>	EPCAM
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates.
<b>Uniprot ID</b>	<a href="#">Q3T0L5</a>
<b>Sensitivity</b>	0.156 ng/ml
<b>Expiration Date</b>	Please enquire.