

## Product Datasheet

### Human IFIT1 ELISA Kit (orb404645)

|                          |  |
|--------------------------|--|
| <b>Description</b>       | Human Interferon-induced protein with tetratricopeptide repeats 1(IFIT1) ELISA kit |
| <b>Reactivity</b>        | Human  |
| <b>Range</b>             | 15.6 pg/mL-1000 pg/mL  |
| <b>Target</b>            | IFIT1  |
| <b>Note</b>              | For research use only  |
| <b>Application notes</b> | Detection Wavelength: 450 nm   |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates, cell lysates                                    |
| <b>Uniprot ID</b>        | <b>P09914</b>  |
| <b>Sensitivity</b>       | 3.9 pg/mL  |
| <b>Expiration Date</b>   | Please enquire.  |