

## Product Datasheet

### Human EB13 ELISA Kit (orb404666)

|                          |   |
|--------------------------|---|
| <b>Description</b>       | Human Interleukin-27 subunit beta(EB13) ELISA kit |
| <b>Reactivity</b>        | Human   |
| <b>Range</b>             | 15.6 pg/mL-1000 pg/mL                             |
| <b>Target</b>            | EB13  |
| <b>Note</b>              | For research use only                             |
| <b>Application notes</b> | Detection Wavelength: 450 nm                      |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates                 |
| <b>Uniprot ID</b>        | <a href="#">Q14213</a>                            |
| <b>Sensitivity</b>       | 3.9 pg/mL   |
| <b>Expiration Date</b>   | Please enquire.                                   |