

Product Datasheet

Human MUSK ELISA Kit (orb404874)

Description	Human Muscle Skeletal Receptor Tyrosine Kinase (MUSK) ELISA Kit
Reactivity	Human
Range	1.56 pg/mL-100 pg/mL
Target	MUSK
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O15146
Sensitivity	0.39 pg/mL
Expiration Date	Please enquire.