

## Product Datasheet

### Human NMU ELISA Kit (orb404922)

<b>Description</b>	Human neuromedin U (NMU) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	0.312 ng/mL-20 ng/mL
<b>Target</b>	NMU
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, cerebrospinal fluid (CSF), tissue homogenates
<b>Uniprot ID</b>	<b>P48645</b>
<b>Sensitivity</b>	0.078 ng/mL
<b>Expiration Date</b>	Please enquire.