

## Product Datasheet

### Human PDCD1LG2 ELISA Kit (orb405069)

|                          |  |
|--------------------------|--|
| <b>Description</b>       | Human Programmed cell death 1 ligand 2(PDCD1LG2) ELISA kit |
| <b>Reactivity</b>        | Human  |
| <b>Range</b>             | 0.45 ng/ml-30 ng/ml.                                       |
| <b>Target</b>            | PDCD1LG2   |
| <b>Note</b>              | For research use only                                      |
| <b>Application notes</b> | Detection Wavelength: 450 nm                               |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates.                         |
| <b>Uniprot ID</b>        | <a href="#">Q9BQ51</a>                                     |
| <b>Sensitivity</b>       | 0.11 ng/ml   |
| <b>Expiration Date</b>   | Please enquire.  |