

Product Datasheet

Human SIGLEC8 ELISA Kit (orb405322)

Description	Human Sialic acid-binding Ig-like lectin 8(SIGLEC8) ELISA kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	SIGLEC8
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9NYZ4
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.