

Product Datasheet

Human SIGLEC9 ELISA Kit (orb405323)

Description	Human Sialic acid-binding Ig-like lectin 9(SIGLEC9) ELISA kit
Reactivity	Human
Range	125 pg/ml-8000 pg/ml
Target	SIGLEC9
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum,plasma,tissue homogenates.
Uniprot ID	Q9Y336
Sensitivity	31.25 pg/ml
Expiration Date	Please enquire.