

## Product Datasheet

### Human SORT1 ELISA Kit (orb405350)

**Description** Human Sortilin(SORT1) ELISA kit

**Reactivity** Human

**Range** 46.88 pg/mL-3000 pg/mL

**Target** SORT1

**Note** For research use only

**Application notes** Detection Wavelength: 450 nm

**Sample Types** serum, plasma, tissue  
homogenates, cell lysates

**Uniprot ID** [Q99523](#)

**Sensitivity** 11.72 pg/mL

**Expiration Date** Please enquire.