

## Product Datasheet

### Human SORL1 ELISA Kit (orb405351)

<b>Description</b>	Human Sortilin-related receptor(SORL1) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	7.8 ng/ml- 500ng/ml
<b>Target</b>	SORL1
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">Q92673</a>
<b>Sensitivity</b>	1.95 ng/ml
<b>Expiration Date</b>	Please enquire.