

Product Datasheet

Human TNFRSF25 ELISA Kit (orb405486)

Description	Human Tumor necrosis factor receptor superfamily member 25(TNFRSF25) ELISA kit
Reactivity	Human
Range	0.24 ng/mL-15 ng/mL
Target	TNFRSF25
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q93038
Sensitivity	0.06 ng/mL
Expiration Date	Please enquire.