

Product Datasheet

Human ADAMTS4 ELISA Kit (orb405624)

Description	human ADAM metalloproteinase with thrombospondin type 1 motif4,ADAMTS4 ELISA Kit
Reactivity	Human
Range	12.5 ng/ml - 800 ng/ml
Target	ADAMTS4
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O75173
Sensitivity	3.125 ng/ml
Expiration Date	Please enquire.