

## Product Datasheet

### Human AP12 ELISA Kit (orb405630)

<b>Description</b>	Human apelin 12,AP12 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	31.25 pg/mL-2000 pg/mL
<b>Target</b>	AP12
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, saliva, urine, tissue homogenates
<b>Sensitivity</b>	7.81 pg/mL
<b>Expiration Date</b>	Please enquire.