

Product Datasheet

Human BAX ELISA Kit (orb405643)

| | |
|--------------------------|---|
| Description | Human Apoptosis regulator BAX ELISA Kit |
| Reactivity | Human |
| Range | 1.25 ng/mL-80 ng/mL |
| Target | BAX |
| Note | For research use only |
| Application notes | Detection Wavelength: 450 nm |
| Sample Types | serum, plasma, tissue homogenates |
| Uniprot ID | Q07812 |
| Sensitivity | 0.312 ng/mL |
| Expiration Date | Please enquire. |